DIGIMAX RF  
TV Listening System

**CLEAR TV AUDIO™**  
*Hear what you see, and love what you hear!*

**Operating instructions**

Page 7

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Página 27

**CLEAR TV AUDIO™**  
*Hear what you see, and love what you hear!*

**CLEAR TV AUDIO™**  
*Hear what you see, and love what you hear!*



Model: \_\_\_\_\_

Owner name: \_\_\_\_\_

Date purchased: \_\_\_\_\_

Distributor name: \_\_\_\_\_

Phone number: \_\_\_\_\_

### Quick overview

These operating instructions allow you to fold out the first 2 pages. These pages show images of the Digimax RF with identification by numbers.

#### Overview headset:

- 1 Volume control
- 2 Loud speaker ear adapters
- 3 Battery compartment
- 4 Status LED
- 5 Pairing button (connect)
- 6 Balance
- 7 Tone adjustment left / right

#### Overview neckloop:

- 8 On/Off switch (On / Off)
- 9 Volume control
- 10 Induction loop
- 11 Breakaway safety connection
- 12 Battery compartment
- 13 Status LED
- 14 Pairing button (connect)
- 15 Balance
- 16 Tone adjustment

#### Overview transmitter:

- 17 Charging cradle
- 18 Charge LED for receiver in charging cradle
- 19 Replacement battery charging bay
- 20 Charge LED for replacement battery charging bay
- 21 Transmitting LED indicator

#### Overview accessories:

- 22 Power supply connection
- 23 USB connection
- 24 AUDIO, MIC / USB selector switch
- 25 Female connector for connection to a sound source (green)
- 26 Female connector for connection to a microphone (black) – optional
- 27 Pairing button (connect)
- 28 Transparent cover
- 29 Digital optical cable (Toslink)
- 30 S/PDIF converter
- 31 L/R audio cable
- 32 Microphone

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<b>Table of contents</b>	<b>Page</b>
<b>Quick overview</b>	<b>6</b>
Overview headphone	6
Overview neckloop	6
Overview transmitter	6
Overview accessories	6
<b>Safety information</b>	<b>8</b>
<b>Package Contents</b>	<b>9</b>
<b>Set up / connections</b>	<b>10</b>
<b>Preferred connections</b>	<b>11</b>
Connecting the transmitter to an RCA female jack	11
Connecting the transmitter to a headphone jack	11
<b>Alternative connections</b>	<b>12</b>
Connecting the TV to the S/PDIF converter	12
Connecting the S/PDIF converter to the transmitter	12
Connecting the transmitter to the microphone (optional)	13
Connecting the transmitter to the USB connector	13
Connecting the transmitter to the power supply	14
<b>How to use</b>	<b>14</b>
Digimax RF with under the chin headset receiver	14
Digimax RF neckloop	15
Safety lock in the neckloop	16
Listening with the 3.5 mm stereo output	16
Pairing of the transmitter and the receiver	17
Balance control	18
Tone control	18
Volume control	19
Transparent cover	19
Range	19
Battery charging	20
<b>Cleaning and care</b>	<b>21</b>
Cleaning the transmitter and the receiver	21
Cleaning the silicone ear adapters	21
Replacing the silicone ear adapters	21
<b>Important information</b>	<b>22</b>
Warranty	22
Environmental information/disposal	23
<b>Troubleshooting</b>	<b>24</b>
<b>Technical Data</b>	<b>26</b>

## Safety information

- Before you begin operating this device, read these instructions carefully and completely.
- Please retain these operating instructions to ensure that they are accessible at any time for other users.
- This device is designed for operation in dry interior spaces, with a room temperature of 50° F to 95° F / 10° C to 35 °C.
- Please handle the device and accessories carefully and install them only in low-dust and clean locations.
- Clean the transmitter and the receiver regularly. For hygienic reasons, the soft silicone ear adapters should be cleaned regularly. Please replace as necessary with the soft silicone ear adapters, offered as special accessories.
- Paints and furniture polishes can degrade the rubber feet of the transmitter and therefore they may cause surface damage to your furniture.
- To power the transmitter, use only the provided power supply.

### Caution:

The receiver volume can be set to a high level. Please comply with the following information to protect your ears or the ears of persons with normal hearing, especially children, from hearing damage due to high volume levels.

1. Set the volume at a low level before you put the receiver on.
2. During use, adjust the volume to a level comfortable to you.
3. Before using the Digimax RF neckloop with a heart pacemaker, please contact your physician or the manufacturer of the pacemaker to ensure that the pacemaker is not subjected to interference.
4. If you use other medical devices, please use the Digimax RF only when adhering to the safety information/safety directives of your physician or the manufacturer.

## Safety information for the lithium polymer battery

The battery used in Digimax RF is a lithium polymer battery, which was produced in accordance with the most current provisions and safety regulations. The battery can be damaged when misused or if not used according to the specifications; in this case Ear Technology Corporation accepts no liability.

### Warning!

The following risks exist in extreme cases:

- Overheating
- Fire
- Explosion
- Smoke and gas

- Do not heat the battery above 158° F / 70° C, e.g., through sun radiation or by casting it into fire.
- Do not subject the battery to moisture.
- Do not short circuit the battery
- Do not disassemble or deform the battery.
- Store device/battery in a cool and dry place at room temperature.
- Use only the provided charger for the battery.
- Recharge the battery in the device regularly by placing the receiver in the transmitter cradle when not in use.
- Charge the battery only at ambient temperatures of 50° F to 104° F / 10° C to 40° C.
- Do not recharge the battery if the device is clearly damaged.
- Please return defective device or devices with worn out batteries to your retailer.


### Package contents

- 1 Digimax RF receiver: under the chin headset receiver or neckloop
- 1 Digimax RF transmitter
- 1 battery AP12A
- 1 screwdriver
- 1 power supply
- 1 audio connection cable with 3.5 mm stereo jack
- 1 audio adapter 3.5 mm to a 1/4" stereo jack
- 1 RCA adapter
- 1 pair of additional silicone ear adapters small (with under the chin headset receiver)
- 1 S/PDIF digital → analog audio converter

Congratulations on the purchase of your Digimax RF system. You have purchased a reliable state-of-the-art device. Please review these operating instructions carefully to ensure that you set up the device correctly and that you are familiar with all functions of the device.

## Set up / connections


Turn off the sound source before connecting the transmitter.

Plug the provided audio connection cable **28** into the green "Audio" female connector **25** of the transmitter.


Several options exist to connect the Digimax RF to a sound source (e.g. TV, computer, etc.).

Please check the output connection options of your sound source.

- |  |   |
|--|---|
|  | – RCA female connector  |
|  | – Headphone connection  |
|  | – External microphone   |
|  | – USB connection (e.g. computer)  |
|  | – S/PDIF digital output—when using digital to analog converter (DAC). If your TV has RCA audio outputs, the included S/PDIF converter is not necessary. |

## Preferred Connections

### Connecting the transmitter to an RCA female jack


Connect the adapter cable with the RCA plugs **29** to the audio connection cable **28**.

Plug the RCA plugs into the RCA connection of your sound source.

- Ensure that the white connection is connected to the white connection of your sound source.
- Ensure that the red connection is connected to the red connection of your sound source.

If your TV has RCA audio outputs, the included S/PDIF converter is not necessary.

### Connecting the transmitter to a headphone jack

Plug the audio connection cable **28** into the headphone jack of your sound source.


You may need to use the 1/4" (6.5 mm) headphone plug adapter (included) to connect to your sound source. Set the volume of your sound source to a comfortable level.

## Alternative Connections

### IMPORTANT: S/PDIF converter


If your TV does not have analog audio outputs, use the S/PDIF converter (included).

### Connecting the TV to the S/PDIF converter

Using the provided Toslink cable **29**, connect your TV (optical audio out) to the S/PDIF converter **30** (optical).

### Connecting the S/PDIF converter to the transmitter


Connect the adapter cable with the RCA plugs **31** to the audio connection cable **28**. Plug the RCA plugs **31** into the RCA connections of the S/PDIF converter **30**.

- Ensure that the white connection is connected to the white connection of your S/PDIF converter.
- Ensure that the red connection is connected to the red connection of the S/PDIF converter.

Plug the 3.5 mm jack of the audio connection cable **28** into the green audio input on the back of the transmitter.

12

### Connecting the transmitter to the microphone (optional)


Use only the A-4967-0 microphone (included).

Plug the microphone **32** into the black "MIC" connection **26** of the transmitter.

Affix the microphone over the loudspeaker of your sound source.

### Connecting the transmitter to a computer using USB connector


The USB cable **33** is not included with your unit. It can be purchased from your electronics retailer. The transmitter has a mini USB connection.

Plug the USB cable into the "USB" connection **23** of the transmitter and connect it to your computer. The switch **24** must be set to "USB".

13

## Connecting the transmitter to the power supply


Connect the provided power supply to the transmitter labeled "DC12V". Plug the power supply into a non-switched electrical outlet.

Charge the receiver with the installed battery for eight (8) hours before using the device for the first time.

## How to use

### Digimax RF with under the chin headset receiver

1. Switch your sound, rce on (e.g. TV).
2. Adjust the volume level control **1** initially to a low volume to protect your hearing. The volume increases when you turn the dial toward "5" and it decreases when you turn the dial in the direction of "0".
3. Put the headset on. The headset switches on automatically when you spread the arms. The green LED **4** indicates that the receiver is switched on.
4. The blue LED bar **21** in the transmitter is illuminated as soon as the transmitter has contact and transmits to the receiver.
5. Set the desired volume **1**.
6. When you remove the headphones they will switch off automatically. We recommend that you put the headphones back into the transmitter cradle when not in use, to ensure that the battery is always charged.


## Digimax RF neckloop

**CAUTION: People using this product who have a heart pacemaker, should contact their physician or the manufacturer of the pacemaker to ensure that the pacemaker is not subject to interference.**

1. Switch your sound source on (e.g. TV), and switch your hearing aids to t-coil setting.
2. Adjust the volume level control **9** initially to a low volume to protect your hearing. The volume increases when you turn the dial toward "5" and it decreases when you turn the dial in the direction of "0".
3. Ensure that the neckloop lock has been closed completely. The neckloop lock is used to comfortably wear the device and to minimize accidents.
4. The receiver will be switched on and off by using the ON/OFF switch **8** on the back. Caution: Always switch the receiver off after use to avoid discharging the battery. The green LED **13** on the receiver indicates that the receiver is switched on.
5. The blue LED bar **21** on the transmitter is illuminated as soon as the transmitter has contact and transmits to the receiver.
6. Set the desired volume.
7. Switch the receiver off after use. We recommend that you put the neckloop receiver back into the transmitter cradle when not in use, to ensure that the battery is always charged.


## Safety lock in the neckloop

(Neckloop version only)


Please ensure that the safety lock is completely closed in the neckloop to guarantee proper functioning.

Please ensure the receiver is worn correctly.


## Listening with the 3.5 mm stereo output

(Neckloop version only)

3.5 mm stereo output connector


earphones

### Important!

The overall operating time of the receiver can be reduced significantly depending on the accessories connected to the output. We recommend the purchase of an additional battery to guarantee longer operating time.

## Pairing of the transmitter and the receiver

If you want to add **additional** receivers (headsets or neckloops) to the transmitter, then the new receivers must first be paired with the transmitter. This is required to ensure that the transmitter and receivers recognize each other and that they can transmit the encrypted data.

This is only required once for each **additional** receiver. A maximum of 4 headsets or neckloops can be paired with a transmitter.

If the receiver does not find a transmitter, then the red LED **4** blinks and no signal will be received.


Please follow the instruction to pair the receiver with the transmitter:

1. Ensure that the transmitter is connected to an electrical outlet and that it is provided with electricity.
2. Put the receiver on and ensure that it is switched on—the green LED should glow. For the Digimax RF neckloop version switch the neckloop receiver on.
3. Press the "Connect" button **27** briefly at the transmitter. The blue LED bar starts flashing.
4. Press the "Connect" button briefly at the headset receiver **5** or neckloop receiver **14** to complete the pairing.
5. Once the pairing has been completed, the blue LED bar on the transmitter and the receiver status LED stop flashing.

### Balance control

(Digimax RF under the chin headset receiver and the 3.5 mm stereo output of the neckloop.)


The balance controls the volume distribution between the right and the left ear with the head set or when the 3.5 mm stereo output is used on the neckloop with earphones, headphones or induction hooks. Adjust the balance so that you have the same volume level in both ears.

The balance adjustment does not work when the neckloop is used with hearing aids.


### Tone control

(Digimax RF under the chin headset receiver and Digimax RF neckloop)


The tone control is used to adjust the frequency response independently in the left and right ears on the headset.

For the neckloop, the (R) tone control adjusts the sound of both the right and left equally.


### Volume control

Set the volume to a comfortable listening level by rotating the volume control from low to high (0-5).

### Transparent cover

A transparent cover (included) may be installed to protect the volume control against unintended adjustments and against contamination by dust or fluids.

Remove both rubber covers and click the transparent cover into the openings.


### Range

Digimax RF uses the newest 2.4GHz technology. In optimal cases the range is up to 145 feet.

The range depends greatly on the building and the environment. The 2.4GHz radio frequency technology also works through walls.


## Battery charging

Once the battery is fully charged, it can be used for up to six (6) hours. The red LED **4/13** starts flashing once the battery is almost empty. The battery should now be charged or exchanged. The operating time is dependent on the listening volume.

After use, put the receiver into the charging cradle of the transmitter, without applying pressure.

Ensure that the receiver is placed correctly. If the receiver is charging, the red LED **18** at the left side of the receiver is illuminated. The LED **18** is illuminated green as soon as the battery in the receiver is charged. For first time operation, it is recommended that the battery be charged for eight (8) hours.

An additional battery can be charged in the transmitter. Additional batteries are available through your Digimax RF distributor.

The additional battery can be plugged into the side charging bay **19**. Please ensure that the sloped side points upward and that you insert it without applying too much force. The LED **20** is illuminated red when the battery is charging. If the LED **20** is green, it is ready for use.

We recommend that you exchange the battery in the receiver with the optional spare battery in the transmitter weekly. This ensures that the batteries are used regularly and that they achieve a maximum service life.

All batteries lose their capacity over time. This is characterized by a shorter usage time.

20


## Cleaning and care

### Cleaning the transmitter and the receiver

If the device becomes dusty or dirty, you can clean it with a soft, clean and damp cloth.

Please disconnect the transmitter from the power supply before cleaning.

### Cleaning the silicone ear adapters

You can remove the ear adapters from the fixture. Clean the ear adapters with mild soapy water. Dry the ear adapters thoroughly and reattach them.


### Replacing the silicone ear adapters

The hypoallergenic silicone ear adapters should be replaced if they become worn or damaged.

New silicone adapters can be ordered from your Digimax RF distributor.

21

## Important information

Product damage caused by improper handling as well as damage due to contamination are excluded from the warranty.

This device is not protected against moisture. Do not place objects that are filled with fluids on or near the device, e.g. vases.

Open fire sources such as burning candles, should not be placed on or near the transmitter.

## Warranty

Refer to the warranty card included in the product package.

## Recycling/disposal

Help us protect the environment. Please take the time to properly dispose of your equipment including the battery. Additional information about recycling of this product can be obtained from your community waste disposal service.

This device adheres to the requirements of the following EU directives:

- 1999/5/EG R&TTE directive
- 2002/95/EG RoHS directive
- 2002/96/EG WEEE directive
- 2004/108/EG EMV directive
- 2006/95/EG low voltage directive

The conformity of the a.m. directives is confirmed by the CE mark on the device. CE conformity declarations are available on the Internet.

Trade Name	Model Number
	Digimax RF
FCC ID: U94SWING1212 This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.	

## Federal Communication Commission Interference Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## Troubleshooting

Symptoms	Cause	Resolution
Receiver charging control lamp is not illuminated, despite the fact that the receiver is inserted in the charging station.	No contact between battery in the receiver and the charging station in transmitter.	Make sure the receiver is in the proper orientation in the transmitter charging cradle.
	No battery inserted into the receiver.	Insert battery in receiver.
Receiver does not work – no sound transmission.	The transmitter is not connected or not connected correctly to the TV set or the stereo system.	Check the cable and the plug connection between the transmitter and the TV or stereo device. Use the operating instructions of the TV or stereo device to check whether the headphone output is correctly connected.
	TV set or stereo system is not switched on.	
	Pairing was lost.	Re-pair transmitter and receiver. (see pg. 17)
	Receiver is not switched on. (Digimax RF neckloop)	Switch the receiver on by using the ON/OFF switch on the back of the neckloop receiver.
	No battery in receiver.	Insert battery in receiver.
	Battery in receiver is not charged.	Charge battery or insert fully charged battery.

Symptoms	Cause	Resolution
The sound transmission is blurred or distorted.	Receiver is too far away from transmitter.	Move receiver closer to the transmitter.
	The battery energy is almost totally consumed.	Charge battery or insert fully charged battery.
The sound transmission is interrupted for a short time and works well afterwards	Receiver is too far away from transmitter. The environment (e.g. walls) distort the reception.	Move the receiver closer to the transmitter.
Battery will only work for a short time.	Battery is not completely charged or defective and therefore no longer operational.	Charge battery or insert fully charged or new battery.
Charge control lamps do not illuminate – despite battery inserted into the charge station of the transmitter.	No power to transmitter.	Check that the power supply is connected to the electrical outlet and the transmitter. Also check that the electrical outlet is working and that no fuse or breaker was triggered.
	No contact between battery and charging station in transmitter.	Insert the batteries correctly and completely into transmitter.

## Technical Data

Modulation process	Stereo (not with neckloop)
Carrier frequency:	2.4 GHz
Tone frequency transmission range:	20 – 20,000 Hz
Distortion factor:	< 1 %
Signal/interference ratio:	> 90 dB
Working range:	32° F – 131° F / 0° C – 55° C

### Under the chin headset receiver:

Operating hours:	approximately 6 h
Maximum volume:	approximately 125 dB
Weight:	approximately 2.04 oz. / 58 g

### Transmitter:

Power consumption:	approximately 4.5 VA
Power supply:	12 VDC
Power supply:	110–120 Volt 50–60 Hz
Weight:	approximately 1.51 oz. / 43 g

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